



Silon-TEX® II is now OEKO-TEX® Standard 100 Certified



Silon-TEX® II is known all over the globe as the premier silicone-lined textile insert material sewn into garments to manage the appearance of hypertrophic scars. Silon-TEX® II is now the only textile insert to pass the rigorous **OEKO-TEX® Standard 100 Certification**. Please visit www.oeko-tex.com for more details.

OEKO-TEX® is a global certification system you can trust. Products with the STANDARD 100 by OEKO-TEX® label have been independently tested for harmful substances according to strict scientific criteria that go beyond global legal regulations

Silon-TEX® II is constructed of medical grade silicone utilizing the patented *Silon-IPN®* technology bonded to a flexible textile material which is certified safe for skin contact. This same technology has been used for over 25 years and now Bio Med Sciences continues to demonstrate itself as a pioneer in manufacturing medical devices with this newest certification distinction.

The benefits of Silon-TEX® II are clear:

- Clinically proven to prevent and reduce scars
- Safe and Effective – over 25 years of successful clinical use in plastic surgery, burn rehabilitation, and dermatology
- Ultra-Thin, Flexible and Breathable – utilizing patented Silon® technology
- Strong and Durable – washable and reusable for effective scar management
- Custom fit to challenging anatomical areas over the entire body as part of a pressure garment
- Proudly made in the USA under strict ISO, FDA and cGMP standards for medical devices

**Silon® Scar Management Products
Proven Safe and Effective for Over 25 Years**



Visit SILON.COM or call 1-800-257-4566.
Outside North America, call +1-610-530-3193 or fax +1-610-530-3194.

Silon is a registered trademark of Bio Med Sciences, Inc. in the United States and/or other countries. This product is protected by one or more U.S. and foreign patents. Additional patents pending. © 2020
Form ID: Silon-TEXIIIEKO-TEX_flyer.indd Rev. 2020-04-30

